INDICAÇÃO:

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INDICATION:

Self-Polymerizable dental acrylic resin, used in Repairs, Additions, Rebases of Denture prostheses and Removable Partial Denture.

CONTRAINDICATIONS:

This medical product must not be used in people with known skin sensitization (contact dermatitis) to polymerizable monomers, as they may cause allergic reactions to hypersensitive people.

WARNINGS/PRECAUTIONS:

- For dental use only.
- Liquid: Flammable and non-toxic material (Flash-point: 10°C/50°F).
- Resin: May form flammable air/dust mixtures.
- As with all acrylic resins, it is recommended to handle the product in ventilated areas, using nitrile gloves and protective goggles; however, handling is not recommended for pregnant women, nursing mothers and people with recognized allergies to this type of product.
- May cause skin irritation in susceptible people.
- Keep away from light and heat. Keep in a dry and cool place.
- Do not dispose of waste into sewers and water courses or in domestic waste.
- Products that are not compatible with the resin/monomer: alcohol, water, chloroform, acetone, petroleum jelly and other solvents.
- The bases of the synthetic teeth must be ready in accordance with the Instructions for Use for the correct support of the prosthesis.
- We do not recommend the use of wax solvents. The residual may cause a lower connection of the teeth with the acrylic base.
- Recommended mechanical polishing for all types of resin.

INSTRUCTIONS FOR USE:

Shake the bottle before using.

- **PREPARATION:** Before starting any job, the containers must be completely clean and free of traces of wax or grease, so that the chemical union between the resins may be perfect.

We recommend a room temperature for mixing and application between 15°C and 30°C.

- **DOSING AND PLASTICIZATION:** Recommended ratio: powder/liquid 2:1 (14g of powder for 7ml of liquid).

Place the necessary amount of liquid in the glass, silicon or Dappen jar, in accordance with the amount of acrylic. Over the liquid, place the resin until the two are mixed together without having excess powder or liquid, if there is, dispose it and homogenize carefully with a spatula, avoiding the incorporation of air bubbles. Observe so that the mixture is liquid, however consistent.

Do not use the resin, when it stops sliding on the spatula (5 – 6 minutes after dosing, at 23°C).

Avoid touching with the fingers on the repaired surface.

- **POLYMERIZATION:** For polymerization, take the resin to a polymerizer, place the necessary amount of water, apply 20 pounds of compressed air and allow it to complete the polymerization, for 20 minutes. Wait for total polymerization to start finishing and polishing. For repairs, additions and rebases, proceed in the conventional way. This resin can be used for orthodontic works. High stability of colors and abrasion resistance.

NOTE - Evoden recommends for this product, solely mechanical polishing.

- **FINISHING AND POLISHING:** Finishing is started with bits and mills, passes to abrasive stones and is concluded with fine textured sandpapers and finishing Rods.

Polishing is performed with soft brush and polishing Rods.

The use of monomer with Crosslink will make polishing and brilliance easier, providing greater durability of the denture.

CHARACTERISTICS/BENEFITS:

- Rigorous selection and raw materials reception tests with successive tests during all manufacture cycles, ensure the results specified for the final product.
- Selection and control of pigments ensures high shade stability.
- Appropriate packaging for use of the supplied measuring utensils, which makes measurement of components easier, separately, and the conservation of the product until the end.
- Ease in handling and application for the execution of the job.

PRESENTATIONS:

POWDER: 50g, 120g, 250g, 500g, 1Kg, 2. 250kg, 5Kg and 50Kg.

LIQUID: 50ml, 120ml, 250ml, 500ml, 1L, 900ml and 5L cans.

SHADES:

Colorless, Pink and Black.

COMPOSITION:**POWDER**

Polymethylmethacrylate
Benzoyl Peroxide
Biocompatible pigments

LIQUID

Methylmethacrylate monomer
DMT
Inhibitor

RECOMMENDATIONS TO PATIENTS:

To increase the durability of the prosthesis, daily hygiene is recommended without the use of aggressive substances, such as acids or alkalis.

SHELF LIFE:

POWDER: 10 years as of the date of manufacture.
LIQUID: 2 years as of the date of manufacture.

ANVISA N°:

POWDER: 82426979003
LIQUID: 82426979006

TECHNICAL RESPONSIBLE:

Jéssica Augusto.
CRQ IV-04270957

ATTENTION: Use our free support channels to obtain the correlation between the purchased product and the corresponding usage instruction version. To obtain, free of charge, these Instructions for Use in printed format, ask the customer service (SAC) on 0800 058 8888.

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