**INDICAÇÃO:**

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**INDICATION:**

Thermopolymerizable acrylic resin developed for prosthetic characterization work according to the patient's gum coloration.

Contraindications: This medical product must not be used in people with known skin sensitization (contact dermatitis) to polymerizable monomers, as they may cause allergic reactions to hypersensitive people.

**WARNINGS/PRECAUTIONS:**

- For dental use only.
- Liquid: Flammable and non-toxic material (Flash-point: 10°C/50°F).
- Resin: May form flammable air/dust mixtures.
- As with all acrylic resins, it is recommended to handle the product in ventilated areas, using nitrile gloves and protective goggles; however, handling is not recommended for pregnant women, nursing mothers and people with recognized allergies to this type of product.
- May cause skin irritation in susceptible people.
- Keep away from light and heat. Keep in a dry and cool place.
- Do not dispose of waste into sewers and water courses or in domestic waste.
- Products that are not compatible with the resin/monomer: alcohol, water, chloroform, acetone, petroleum jelly and other solvents.
- We do not recommend the use of wax solvents. The residual may cause a lower connection of the teeth with the acrylic base.
- Recommended mechanical polishing for all types of resin.

**INSTRUCIONES FOR USE:**

To use the dosing nozzle, remove the seal and cap completely.

Before starting any job, the containers must be completely clean and free of traces of wax or grease.

The universal characterization system for dental prosthetics is very versatile and can be used for various polymerization cycles. Here are 3 examples of polymerization to achieve a residual monomer content of less than 1%:

**CONVENTIONAL POLYMERIZATION:** Place the muffle in a pot with 3 liters of cold water, turn on the stove with low flame, allowing it to reach a temperature of 70°C (158°F). Maintain this temperature for 30 min. Increase the flame to reach 100°C (212°F) keeping it boiling for 90 minutes and then, put out the flame. Allow the muffle to cool down inside the water to 40°C (104°F) for approximately 20 minutes, before starting the deflasking.

- **THERMO-PNEUMATIC POLYMERIZATION:** Use a polymerizer with pressure and temperature indicators, insert the flask covering it with cold water. Close the lid and apply the following instructions:
- Inject 60 pounds/cm<sup>2</sup> of compressed air;
- Connect the polymerizer until it reaches a temperature of 100° C;

- Turn off the polymerizer waiting until the temperature drops to 40° C;
- Reconnect until the temperature reaches 100° C again; • Turn off the polymerizer and allow it to cool naturally.
- **THE MICROWAVE POLYMERIZATION** technique is indicated for pieces of 10mm maximum thickness, using Termoden liquid with Crosslink, the use of silicone barrier is recommended.

Microwave polymerization cycles may vary depending on the power of each device. Check how to program, increase and/or decrease the power of your device from instructions manufacturer.

How to know the power of a microwave device? Check the inside edge of the microwave rack for a label that specifies wattage. If you don't know, look for the serial number on the

back of the device that can show power data.

It is important to control and preserve the muffle moisture.

In cases of work that will be carried out overnight, keep the muffle hydrated.

To obtain complete polymerization, it is recommended to follow the parameters defined in the table below:

Microwave 500 Watts	
Initial Stage	20 min. - Power 20/30%
Final Stage	+5 min. - Power 80/100%
When there are 2 muffles + 10% power in the final cycle	

Microwave 800 Watts or 1100 Watts	
Initial Stage	20 min. - Power 10/20%
Final Stage	+5min. - Power 50/60%
When there are 2 muffles + 10% power in the final cycle	

Microwave 1200 Watts or 1400 Watts	
Initial Stage	20 min. - Power 10%
Final Stage	+5 min. - Potência 30/40%
When there are 2 muffles + 10% power in the final cycle	

Note: If the microwave does not have 0% power, use the lowest power.

**CHARACTERISTICS / BENEFITS:**

- The method allows the customization of the prosthesis according to the client's gum coloration.
- It supports the most varied polymerization methods for professionally known thermosetting resins.
- Low water absorption, low residual monomer, excellent mechanical properties, high impact strength, unequaled gloss, absence of porosities, easy flow and compaction.

**PRESENTATIONS:**

Powder: Bottle 50g

Kit with:

- 9 bottles of resin powder of 25g each;
- 2 bottles of veins.

**SHADES:**

Light Pink, Medium Pink, Red Pink, Purple Pink, Brown Pink, Black Pink, Intensive Red I, Intensive Red II, and Intensive Purple.

**COMPOSITION:**

POWDER: Polymethylmethacrylate, Benzoyl Peroxide, Biocompatible pigments

VEINS: Rayon yarns and pigments

**SHELF LIFE:**

10 years as of the date of manufacture.

**ANVISA N°:**

82426979001

**TECHNICAL RESPONSIBLE:**

Jéssica Augusto.  
CRQ IV-04270957

**ATTENTION:** Use our free support channels to obtain the correlation between the purchased product and the corresponding usage instruction version. To obtain, free of charge, these Instructions for Use in printed format, ask the customer service (SAC) on 0800 058 8888.

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