INDICAÇÃO:

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INDICATION:

It is used in the machining technique for making ceramic casting and injection patterns (Burnout) - machined by CAM System from various manufacturers.

CONTRAINDICATIONS:

Contraindicated for people allergic to any of the components.

INSTRUCTION FOR USE:

Attach the block to the machine following the equipment manufacturer's instructions.

CHARACTERISTIC AND BENEFITS:

- Excellent grindability
- Excellent translucency
- Excellent stability
- Excellent strength
- Burns without residue

MECHANICAL PROPERTIES

Modulus of Elasticity	3.050 MPa
Bending Strength	96,55 MPa
Toughness (Charpy): (Izod):	1,52 KJ / m ² 24,70 J / m
Tensile Strength	72,19 MPa
Vickers Hardness HV (0,05/10)	17,33
Water Sorption ISO 10477	Max 32 µg/mm ³
Solubility in water	Max 1,6 µg/mm ³
Residual monomer:	Max 2,2 %

MILLING PATTERNS FOR BURNOUT

Please use the suggested furnace cycle/program shown below:

STEP HEATING	HEAT RATE	X1	X3	X6	X9
Room temperature to 260°C (500°F)	3°C (37.4°F)/min				
Holding time at 260°C (500°F)*		40 min.	60 min.	70 min.	90 min.
260°C to 580°C (500°F to 1076°F)	6°C (42.8°F)/min				
Holding time at 580°C (1076°F)*		20 min.	30 min.	40 min.	50 min.
580°C (1076°F) to end temperature**	9°C (48.2°F)/min				
Holding time at end temperature*		30 min.	40 min.	50 min.	60 min.

Important Notes:

* If more rings are inserted together in a furnace, each holding time should be extended by 10 minutes.

** Use the end temperature indicated by the investment and restoration material being used, as metal (casting), ceramic pressing, etc.

MODELS:

EBW, EBZ e EBA.

COLORS:

Red and Colorless

HEIGHTS:

According to the table

Models/ Heights (mm)	13	14	16	20	22	25
EBA	X	-	-	X	-	-
EBW	-	X	X	X	X	-
EBZ	-	X	X	X	X	-

COMPOSITION:

Polymethyl methacrylate, EDMA, Fluorescent and Biocompatible Pigments.

EXPIRATION DATE:

10 years after the manufacturing date.

PACKAGING:

1 piece.

STORAGE:

keep in a dry place and protected from sunlight.

ANVISA N°:

82426979008

TECHNICAL RESPONSIBLE:

Jéssica Augusto.
CRQ IV-04270957

ATTENTION: Use our free service channels to obtain the correlation between the supplied product and the version of the corresponding Instructions for Use. To obtain, free of charge, these Instructions for Use in printed form, please ask Customer Service at 0800 058 8888

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