INDICAÇÃO:

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INDICATION:

Cobalt chromium alloy for removable partial dentures.

CONTRAINDICATIONS:

Sensitive patients may have hypersensitivity reactions to one or several of the metals contained in the alloy.

WARNINGS/PRECAUTIONS:

- The product is intended exclusively for dental applications.
- Keep the product out of humidity and out of reach of children. Sensitive patients may have hypersensitivity reactions to one or more of the metals contained in the alloy.
- In case of skin contact with casting residues or vapors, wash with abundant water and neutral soap.
- The reuse of cast materials is not recommended.

INSTRUCTIONS FOR USE.**INCLUSION:**

Evolay Alloy is compatible with all coatings sold in the market.

Evolay indicates the use of Evolay Coating, which requires a preheating of 850°C.

FOUNDING:

Evolay alloy can be cast with all casting equipment.

INDUCTION FOUNDRY:

When using induction casting equipment, do not use graphite or carbon crucibles. Pour before the oxide layer is broken.

- After casting, the ring must cool at room temperature for at least 20 minutes, do not immerse in water.
- Blast the object with aluminum oxide or similar materials.
- To avoid destroying the passivation layer on the surface, the alloy must not be treated with pickling agents.

FINISHING:

Evolay Cobalt-Chromium alloy for removable partial dentures can be machined with any usual instrument for finishing cast prostheses.

WELDING:

Weld the bonding spaces of the Evolay Cobalt-Chromium alloy with precious

alloys. The solder used must be selected according to the precious metal to be soldered.

REUSE:

Evolay does not recommend the use of leftovers as it may contaminate the entire process.

CHARACTERISTICS/BENEFITS:

- The rigorous selection of raw materials, receiving tests combined with successive tests during all manufacturing cycles, guarantee the specified results in the final product.
- Evolay alloy has excellent strength and toughness, being wear and corrosion resistant.
- It is compatible with all coatings commercialized in the market.
- It can be worked with any usual instrument for finishing cast prosthesis.

PRESENTATIONS:

Cylinders of approximately 8.3g each and are stored in 50g packages.

COMPOSITION:

Cobalt 60%, Chromium 30%, Molybdenum 5%, Tungsten <1%, Silicon <1%, Manganese <1%, Nickel <1%.

STORAGE:

Keep in a dry and airy place, protected from light, humidity and contaminants.

VALIDITY:

Undetermined expiration date

ANVISA N°:

82426979009

TECHNICAL OFFICER:

Jéssica Augusto
CRQ IV-04270957

ATTENTION: Use our free support channels to obtain the correlation between the purchased product and the corresponding usage instruction version. To obtain, free of charge, these Instructions for Use in printed format, ask the customer service (SAC) on 0800 058 8888.

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