**INDICAÇÃO:**

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CRQ IV-04270957

**RECOMMENDATION:**

Evotray is a light-cured impression mold material in the form of prefabricated base plates used to create custom rigid impression trays and bite blocks.

- Custom impression tray.
- Bite wafer for try-in framework.
- Key for supporting pin registration.
- Bite -block.
- Splinting outer coping or telescoping crowns.
- Baseplate for setting up dentures.

**CONTRAINDICATIONS:**

The use of this medical product is contraindicated in people with known skin sensitization (contact dermatitis) to polymerizable monomers, as they may cause allergic reactions in hypersensitive people.

**USE INDICATIONS:****PREPARATION**

Mark the outline of the future tray on the model.

Mask the undercuts on the model.

To create spacer for the impression material, coat the entire model by means of a sheet of wax.

If the plaster is very dry, sprinkle the model slightly with water or apply a alginate insulation.

**APPLICATION**

- Take the baseplate from the box and close it.
- Once the model has been coated, adapt the baseplates and cut to shape.
- Use the remaining material to form a handle and attach it to the tray.
- Place the model and tray in a light-curing unit and cure without vacuum.
- Remove the tray from model and polymerize at the inside.
- Use alcohol to remove the smear, tacky layer.

**WORKING TIME**

- Working in daylight: 20 minutes
- Polymerization in an ultraviolet light-curing unit
- The top of the tray, on the model: 3-5 minutes.
- The inside of the tray, without model: 3-5 minutes.

**TRIMMING**

Always use a dust extractor when trimming the tray!

Trim the custom tray with tungsten carbide burs and smooth it with sandpaper.

Pumice slurry can be used for polishing.

**WARNINGS AND PRECAUTIONS:**

When there is a proven allergy to any of the components of the formula, the product should not be used.

Use the product in a ventilated place, wearing gloves, mask and protective glasses.

Handling this product is not recommended for pregnant and breastfeeding women.

Keep the packaging closed.

Do not reuse packaging.

Correctly follow the instructions for use, in order to maintain the product properties and performance defined by the manufacturer.

Follow the usual procedures for microbiological control during handling.

Do not dispose of waste in domestic sewers.

For dental use only.

Store away from light and heat.

**TECHNICAL DATA:**

Flexural strength : 120N/mm<sup>2</sup>

**PACKAGING:**

Packaging with 10 units- colorless

Packaging with 50 units- colorless

**INGREDIENTS:**

Tripropylene glycol diacrylate; Phosphine oxide

**SHELF-LIFE**

3 years from the date of manufacture.

ANVISA N°: 82426970003

**TECHNICAL RESPONSIBLE:**

Jéssica Augusto

CRQ IV-04270957

**ATTENTION:** Use our free support channels to obtain the correlation between the purchased product and the corresponding usage instruction version. To obtain, free of charge, these Instructions for Use in printed format, ask the customer service (SAC) on 0800 058 8888.

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