

UPPER AND LOWER
POSTERIOR

10°		10/30 (L5) 30 (L6) 7,7 (L7) 32 (L8) 8,25
10°		10/32 (L5) 32 (L6) 8,5 (L7) 34,5 (L8) 9,25
10°		10/34 (L5) 34 (L6) 9,1 (L7) 36 (L8) 9,5
20°		29M (L5) 29 (L6) 8,25 (L7) 32,5 (L8) 9,1
20°		31M (L5) 31 (L6) 8,75 (L7) 34,5 (L8) 9,7
20°		33M (L5) 33,2 (L6) 9,25 (L7) 37 (L8) 10,2



0800 058 8888
sac@evoden.com.br
evoden.com.br

ACRYLIC RESIN TEETH FOR DENTAL PROSTHESES
Follow the usual procedures for choice, waxing-up, articulation, polymerization, infection-control and maintenance of the prosthesis, as well as for individual protection.
Indications for the patient: To increase the teeth durability it is highly recommended to avoid contact with alkalis, acids or other substances that can damage acrylic materials during the daily hygiene.
Composition: Polymethylmethacrylate - Dimethacrylate - Fluorescent - Biocompatible Pigments.
Storage: Keep protected from light and heat.
Packing Units: 8 anterior sets or 5 posterior sets.

CERTIFIED COMPANY



EVODEN INDÚSTRIA, COMÉRCIO, IMPORTAÇÃO E EXPORTAÇÃO DE PRODUTOS ODONTOLÓGICOS LTDA.
Rua Genésio de Bem, 1634, Distrito Industrial.
CEP 13.633-428. Pirassununga, São Paulo, Brasil.
CNPJ: 44.136.221/0001-67. Indústria Brasileira.
Made in Brazil. Hecho en Brasil.
Fone + 55 (19) 3565-8888.
Anvisa nº: 82426979002
Responsável Técnico:
Jéssica Augusto. CRQ IV-04270957.

Catálogo Dental - 50789 - Rev. 00



EVODEN



ACRYLIC TEETH **2** LAYERS
MOLD CHART

UPPER ANTERIORS

2D
(L1) 41,0
(L2) 7,4
(H1) 9,2
(3) 10,4

3D
(L1) 43,7
(L2) 8,1
(H1) 9,0
(3) 10,6

2N
(L1) 40,8
(L2) 7,5
(H1) 9,8
(3) 11,7

3N
(L1) 42,6
(L2) 8
(H1) 10
(3) 11,7


3M
(L1) 38,4
(L2) 7,1
(H1) 9,1
(3) 10,4

3P
(L1) 47,6
(L2) 8,4
(H1) 10,3
(3) 11,8

A25
(L1) 43,2
(L2) 8,4
(H1) 9,2
(3) 11,1

A26
(L1) 45,8
(L2) 8,7
(H1) 9,4
(3) 11,1

263
(L1) 42,1
(L2) 7,6
(H1) 9,7
(3) 11,7



TOOTH SHAPE
ACCORDING TO BIOTYPE



TAPERED
SLENDER



SQUARE
ATHLETIC



OVAL
ROBUST



UPPER ANTERIORS

264
(L1) 44,1
(L2) 8
(H1) 10,3
(3) 12,3

266
(L1) 49
(L2) 9,1
(H1) 11,7
(3) 14,7

38
(L1) 49,4
(L2) 9,2
(H1) 11
(3) 12,8

LOWER ANTERIORS

2D
(L3) 31,3
(L4) 4,9
(H2) 9,0
(3) 10,7

3D
(L3) 32,8
(L4) 5,1
(H2) 8,7
(3) 10,1

2N
(L3) 32,4
(L4) 5,0
(H2) 8,8
(3) 11,4

3N
(L3) 33,5
(L4) 4,9
(H2) 9
(3) 10,9

3M
(L3) 31,1
(L4) 4,7
(H2) 8,4
(3) 9,8

3P
(L3) 37
(L4) 5,6
(H2) 9,7
(3) 11,4

A25
(L3) 34,6
(L4) 5,4
(H2) 9,4
(3) 11,1

A26
(L3) 35,8
(L4) 5,4
(H2) 8,4
(3) 10,5

ARTICULATION CHART

ANTERIORS	0°	10°	20°
2D	0/32	10/32	29M
3D	0/32	10/32	31M
2N	0/32	10/30	29M
3N	0/32	10/32	31M
3M	0/32	10/30	29M
3P	0/32	10/34	33M
A25	0/32	10/32	31M
A26	0/32	10/32	31M
263	0/30	10/30	29M
264	0/32	10/32	31M
266	0/32	10/34	33M
38	0/32	10/34	33M

SHADE GUIDE

BL1 - BL4 - A1 - A2 - A3 - A3,5
B1 - B3 - C2 - C3 - D3

DIFFERENTIALS



PRODUCTION
TECHNOLOGY BY
CAD/CAM SYSTEM



SYMMETRICAL
PAIRS



ORGANIC ADDITIVE - ODT
*Increases wear and impact
resistance.*



ESTHETIC
*Textured surface, evidenced
mamelons and anatomy that
bring an aesthetic and natural
result to the teeth.*




CONSTANT QUALITY
STANDARD

263
(L3) 32,4
(L4) 5,0
(H2) 8,8
(3) 11,4


264
(L3) 34,6
(L4) 5,4
(H2) 9,4
(3) 11,1

266
(L3) 37,6
(L4) 5,6
(H2) 11,4
(3) 12,9

38
(L3) 37,6
(L4) 5,6
(H2) 11,4
(3) 12,9



UPPER AND LOWER
POSTERIOR



POSTERIOR ANGULATIONS



0°

0/30
(L5) 30
(L6) 8,8
(L7) 32
(L8) 8,8

0°

0/32
(L5) 32
(L6) 9,4
(L7) 34
(L8) 9,4